

LASER TRIANGULATION OF THE FEMORAL HEAD
FOR TOTAL KNEE ARTHROPLASTY ALIGNMENT
INSTRUMENTS AND SURGICAL METHOD

BACKGROUND OF THE INVENTION

5 The present invention relates to laser triangulation of the femoral head for total knee arthroplasty alignment instruments and surgical method. Orthopedic surgeons have been struggling with the alignment of total knee arthroplasties since their inception in the early 1970s. Basically, what is necessary is a 5-7 degree
10 angular resection of the distal femoral condyles as related to the mechanical axis of the femur and a perpendicular resection of the proximal tibia as related to its central axis. The **mechanical axis** is defined as a line extending from the center of the femoral head through the center of the knee to the center of the ankle. Early
15 on, resections of the distal femur and proximal tibia were made by visually trying to match the existing anatomy by eye. Alignment varied considerably depending on the skill of the operating surgeon.

20 In the early 1980s, precision jigs were introduced that aligned the resections to the mechanical axis of the respective bones. Initially, these were **extramedullary** jigs, where the alignment was done outside of the bone by extended rods with perpendicular cutting heads. Later, the **intramedullary** rod method

became popular. This technique required the rod to extend through the intramedullary canal of the patient's femur and tibia.

Many sophisticated instrument systems have emerged over the last 25 years to make these difficult resections more reliable and reproducible for the average orthopedist. Approximately 250,000 total knee arthroplasties are done in the USA each year. After 35 years, there is a wealth of knowledge of the failure modes of this procedure. The most important parameter is accurate alignment of the components. It has been proven that only 4.5 degrees of misalignment causes the components to only load one side of the knee joint leading to rapid failure of the implant. The literature strongly supports the conclusion that the closer the surgeons approach neutral alignment, as defined as equal weight bearing on both sides of the joint, the more successful the implant system will be with longevity, so far, up to 30 years in some reported series. Misaligned total knee arthroplasties tend to get worse with time because the abnormal weight distribution accelerates the wear on the overloaded side leading to rapid failure within a few years in the case of the grossly malaligned.

Although total knee alignment is made up by two resections, one on the distal femur and one on the proximal tibia, the proximal tibia resection is relatively easy because the landmarks are visible to the surgeon. Consequently, this discussion will be

limited to the **distal femoral cut**, which is blind to the surgeon because the femoral head location cannot easily be determined.

There are basically two standard alignment methods for total knee surgeries. The first, as explained above, is called the **extramedullary method** meaning the alignment is accomplished without inserting a rod into the intramedullary canal of bone. This method requires the femoral head to be located either by x-ray or the surgeon's educated guess. Then an extended rod with a perpendicular resection guide is used to make the cuts on the distal femoral condyles by holding one end of the rod at the center of the femoral head and the perpendicular resection guide at the level of resection at the center of the knee.

Two patented systems known to Applicant disclose a variation of this **extramedullary method**. The first is a patent by Dance et al. (#5,690,638) called "A Method and Apparatus for the Alignment of a Femoral Knee Prosthesis." This method uses a pivotable resection guide at the distal condyles that allows a distal retraction force in a freely suspended knee that locates the center of rotation of the femoral head. The resection guide is then locked into position and the cuts made on the distal femur.

There are three problems with this system. The first is that it requires a rigid fixation system of the resection head to the distal femur because a substantial force is necessary to rotate the cutting head to the center of rotation. Their described attachment

system requires extensive resection of soft tissues from the distal femur. The second problem is that the system requires the patient's entire leg be suspended freely so no external restrictions will affect the subsequent traction procedure. The elaborate suspension system is expensive and time consuming to set up. The third problem is that the surgeon has no way to check the accuracy of the traction procedure and has to proceed on faith that the system is accurate.

Another variation of the extramedullary method is disclosed in a patent disclosing an invention for which the applicant herein, Thomas D. Petersen, is a co-inventor, (#5,606,590), titled "A Surgical Laser Beam-Based Alignment System and Method." This patent is an improvement over the extramedullary alignment assist device disclosed in U.S. Pat. No.4,524,766, "A Surgical Knee Alignment Method and System." In this prior patent, the use of radiographically opaque L arm located over the patient's femoral head was described.

Applicant's co-invention (#5,606,590) utilizes a laser and a sophisticated X-ray cassette with a radiograph scale above and below the patient's hip to enable the surgeon to minimize parallax error. If the X-ray beam is not exactly perpendicular to the patient's hip, there is parallax error showing up on the X-ray as to the location of the femoral head. A laser is then set up in line with the center of the patient's femoral head and used to fix

the distal femoral resector perpendicular to this longitudinal axis.

The problems with this system include the elaborate X-ray system that requires proper positioning of the patient on the custom X-ray cassette and the time it takes to set up, take and interpret the X-ray. Then the laser needs to be adjusted to the center of the femoral head and calibrated to the longitudinal axis. Although everything has been done to eliminate parallax error, the patient's femur is only roughly equidistant between the radiographic scales so there is still residual inherent error present in the system, albeit only more than one degree in 10 percent of patients.

The other, more commonly used, method of alignment for total knee arthroplasties, as explained above, is the **intramedullary method**. This method uses a rod placed through a drill hole in the notch between the distal femoral condyles into the intramedullary canal of the long axis of the femur. The distal end of the rod has a protractor adjustment that allows the surgeon to dial in an angulated resection of the distal femoral condyles, typically 5-7 degrees laterally off the central axis of the rod.

There are three major problems with intramedullary alignment systems. The first is an inherent spatial mechanical error of up to 2-3 degrees in both sagittal and axial alignment because the intramedullary canal has a cross-section about the size of a

nickel. Consequently, the rod can be angulated within that space in both planes, worse case, 2-3 degrees, especially in larger femurs.

The second major problem is, not infrequently, the femur is deformed so that it cannot provide a reliable guide. Congenital bowing seen in dwarfism such as Morquio's disease and bone dysplasias are examples where this system should not be used. Previous trauma with resultant angulations of the canal is another example where the intramedullary system should not be used because of significant error.

Lastly, there are well-documented medical complications when the **intramedullary** system is used that are not present with the extramedullary approach. The medical literature is full of articles that document the increased risk from fat embolism and increased blood loss with the intramedullary method. Usually, there is 50 percent more blood loss when the intramedullary canal is used for alignment. Patients require more blood transfusions with this method and accept the risks that go with blood transfusions.

Fat embolism can be a serious problem. The intramedullary fat, dislodged by the drilling and placement of the rod, gets into the patient's blood stream and in most cases causes lethargy and sensorial changes that inhibit the patient's post-operation rehabilitation and delay hospital discharge at least one day.

With these problems in mind, the present invention was developed to improve alignment for total knee arthroplasties and subsequent increased longevity of the implant system while minimizing the intrusiveness of the instruments.

SUMMARY OF THE INVENTION

The present invention relates to laser triangulation of the femoral head for total knee arthroplasty alignment instruments and surgical method. The present invention consists of a newly developed extramedullary method to align the Distal Femoral Resection for Total Knee arthroplasty without violating the intramedullary canal. This method preferably uses a miniature high tech laser to triangulate the center of the femoral head without the need for an inter-operative X-ray.

By avoiding IM rods for alignment, considerable morbidity will be spared. In studies by Applicant using Applicant's Laserlign I system, covered by U.S. Patent No. 5,606,590, an extramedullary alignment system in over 500 cases (several surgeons over the past 5 years) has shown marked improvement in the following: 1) 50% lower blood loss rarely requiring transfusion; 2) marked improvement in the patient's sensorium and lethargy post surgery due to subacute fat embolism, well documented in several studies, allowing the patients to be discharged from the hospital, on the

average, one day earlier; 3) improved alignment greater than 90% within 1.5 degrees of optimal alignment.

Surgeons who have used the Laserlign I system appreciate the Laser Ruler that can be used throughout the case to confirm the Mechanical Axis allowing for inter-operative adjustments.

The present invention, known as the Laserlign II, is actually more accurate than the Laserlign I system that required an X-ray to determine the center of the femoral head. It is intuitive and based on sound mechanical principles that leave no doubt the analysis is correct. The laser gives the surgeon visual confirmation compared to other extramedullary systems that are based on faith that the mechanical principle is correct.

The Laserlign II system can be used with any commercial Total Knee System. The Laser keys up on the Distal Femoral Resector, therefore this vital instrument is included in the Laserlign II system. Since the Laser/Resector must be as flat as possible to the anterior femoral cortex, several other instruments that prepare the distal femur for the Laser/Resector are included as well. The other commercial systems pick up after the preliminary cut on the anterior femoral condyles and the Laserlign II resection of the distal femoral condyles.

The procedure includes a Laser that is set up over the femoral head on a platform that creates a Laser Ruler, i.e., the Mechanical Axis, for the surgeons to refer to throughout the case to make sure

their alignment is correct. It is extremely gratifying to be able to confirm on the operating table that the alignment is correct.

The present invention utilizes a small diode laser positioned at the center of the patient's knee that is adjustable to the longitudinal axis of the femur to facilitate triangulation of the center of the patient's femoral head.

The present invention includes the following interrelated objects, aspects and features:

(1) The present invention includes a removable diode laser that attaches to the distal femoral resector guide at the knee. A V-Frame positioning instrument is provided that has a central pivot allowing adjustment of the longitudinal axis of extension of the laser beam about this central pivot.

(2) The distal femoral resector guide is attached by a pivot post on the V-Frame, and a removable retractor shield is attached thereto.

(3) A clamping mechanism is provided to facilitate clamping of a laser target in a suspended location above the surgical table. In this regard, a clamping device clamps a vertical attachment arm to the table and a horizontal bar supports the laser target.

(4) The method of performing surgery in accordance with the teachings of the present invention is disclosed in detail. After the laser target has been positioned over the hip, an estimate is made of the location of the femoral head. Thereafter, the distal femur is prepared including preparation of a centralized

intercondylar notch and a flat plane on the anterior femoral cortex using round and flat rasps, respectively. Thereafter, the top portions of the anterior femoral condyles are removed using a resector guide for that purpose.

5 (5) The inventive V-Frame is attached to the distal femur and a femoral resector guide is attached to the V-Frame. With this instrumentation in place, the next step is to mount the laser beam generator in a tapered slot on the distal femoral resector guide.

10 (6) Surgical techniques to be described in greater detail hereinafter facilitate location of the center of the femoral head using the laser beam on the target that had previously been located in an estimated location with regard to the femoral head. When these adjustments are completed, the bulls-eye of the laser target is now positioned directly over the center of the femoral head and
15 the laser beam extends precisely down the mechanical axis of the leg.

 (7) With these steps being accomplished, accurate resection of the distal femoral condyles may be performed to facilitate completion of total knee surgery.

20 As such, it is a first object of the present invention to provide laser triangulation of the femoral head for total knee arthroplasty alignment instruments and surgical method.

It is a further object of the present invention to provide such an apparatus and method in which a target is positioned over the hip of the patient.

It is a still further object of the present invention to provide such a device, which allows interconnection with a distal femoral resector guide.

It is a yet further object of the present invention to provide such a device in which a laser beam is used to dramatically enhance the accuracy of the resections that must be carried out during the performance of total knee surgery.

These and other objects, aspects and features of the present invention will be better understood from the following detailed description of the preferred embodiment when read in conjunction with the appended drawing figures.

BRIEF DESCRIPTION OF THE DRAWINGS

Figure 1 shows a side perspective view of a surgical table with a patient thereon and a target positioner mounted on the surgical table.

Figure 2 shows a further side perspective view rotated from the view of Figure 1, partially exploded, and showing the laser target mounted on the attachment arm.

Figure 3 shows a further side perspective view depicting the manners of adjustment of the laser target with respect to the hip of the patient.

Figure 4 shows a perspective view of a distal femur showing the position of a rasp used to create a centralized intercondylar notch.

Figure 5 shows three separate views of the distal femur showing the use of a flat rasp to create a flat plane on the anterior femoral cortex of the distal femur.

Figure 6 shows the use of an anterior condyle resector guide to resect a planar surface on the top of the distal femur by resecting portions of the anterior femoral condyles.

Figure 7 shows a front perspective view of the distal femur showing the mounting of a V-Frame thereon.

Figure 8 shows a side perspective view of the distal femur with a distal femoral resector guide mounted on the V-Frame.

Figure 9a shows an exploded side view of a laser device.

Figure 9b shows a cut away-section of the laser device of Figure 9a.

Figure 10 shows a side perspective view of the laser device mounted on the distal femur via the distal femoral resector guide.

Figure 11 shows use of the laser to align a planar laser beam with the target mounted over the hip of the patient.

Figures 12 and 13 show a front perspective view of the Target illustrating the geometry of alignment of the knee laser beam to the center of the femoral head.

Figure 14 is a blown-up of the front perspective view that illustrates correction of the knee laser beam to the center of the femoral head.

Figure 15 shows a front perspective view of an adjustment knob used to align the knee laser beam.

Figure 16 shows a front perspective view of the laser target when the target is correctly placed over the center of rotation of the femoral head.

Figure 17 shows a side perspective view of the laser beam properly aligned along the mechanical axis of the patient's leg.

SPECIFIC DESCRIPTION OF THE PREFERRED EMBODIMENT

Reference is first made to Figures 1-3, which depict the initial positioning of a target over the hip of the patient. The patient is designated by the reference numeral 1 and has a leg 2 including an upper leg 3, a lower leg 4, and a knee 5, as well as a hip 6 (Figure 3). A surgical table 10 has a surface 11 on which the patient 1 is positioned. A non-sterile target positioner 13 includes a non-sterile clamp 15 that is suitably clamped to the table 10, a vertical portion 17, and a horizontal portion 19 from

which extends a perpendicular portion 21 designed to engage the patient's thigh when the upper leg 3 is flexed 90 degrees with respect to the patient's body (Figure 1) and the surgical table, to position the patient on the operating table. Once the patient is properly positioned, the target positioner is removed from the clamp 15, whereupon the laser target system 25 is mounted on the clamp 15 (Figures 2-3).

As particularly seen in Figures 2 and 3, the laser target system 25 includes a bracket consisting of an attachment arm 18 sized to be received by the clamp 15 and a horizontal slide 27 on which the target portion 29 with a complementary recess 31 is mounted. As seen in Figures 2 and 3, the interaction between the horizontal slide 27 and the recess 31 allows sliding movements in the direction of arrow 33 (Figures 2 and 3). Additionally, the target portion 35 is mounted on a further slide mechanism 37 allowing reciprocation of the target 35 in the direction of the arrow 39 (Figure 3).

The vertical support 18 has the horizontal slide 27 mounted thereon using a mount 41 that has an opening through which the vertical support 18 extends and a screw clamp 43 that allows fixing the vertical position of the mount 41 in any desired vertical position. Thus, through use of the adjusting means described above, the target portion 29 may be adjusted in position in three

degrees of freedom, all the movements being perpendicular to each other, taking into account all the adjustments described above, to allow the target portion 29 to be positioned over the hip 6 of the patient 1 so that, during the course of practice of the inventive surgical procedure, a laser beam 107 (Figure 11) can be aligned over the target 35 while also extending precisely along the mechanical axis of the patient's leg. As seen in Figures 2-3, the target portion 29 has a longitudinal arm comprising a bulls-eye 36 on its top surface and a front surface comprising a vertically depending surface 38 having a plurality of vertical lines 42 (Figures 3 and 11) provided to facilitate alignment of the laser.

With reference now to Figures 4-6, an explanation will be made of the steps that are taken in practicing the inventive method to prepare the distal femur for attachment of a V-Frame and a distal femoral resector guide.

As shown in Figure 4, a rasp 45 having a round cross-section is employed to centralize the intercondylar notch 47 located between the condyles 44 and 46 forming the filed notch 48 (Figure 5). Thereafter, with reference to Figure 5, a flat rasp 49 is used to create a flat plane 52 on the anterior femoral cortex 51 of the distal femur. Thereafter, with reference to Figure 6, an anterior condyle resector guide 53 is placed on the planar surface 52 formed

by the rasp 49, with the guide 53 having a flat guiding surface 55. A surgical saw 57 is employed to resect the top portion of the anterior femoral condyles as shown in the right-hand view of Figure 6, to form a planar surface 59, and a planar surface 61, with the surfaces 59 and 61 being 1/8" higher than the surface 52 but parallel with the surface 52 formed by the rasp 49. This is seen in Figure 6.

With reference to Figure 7, a V-Frame 65 includes a downwardly depending pair of angled brackets 67 and 69 that are attached to the distal end of the femur using cortical screws 71. The V-Frame 65 includes angled wings 75 and 77, and a removable tongue 79 extends proximally toward the hip lying on the plane formed by the surface 52 (see Figure 5).

The V-Frame 65 also includes an upwardly extending post 81 at the apex between the wings 75 and 77. Wing 75 has engraved rotational indicia 78 on its top surface 78. Extending downward is the notch guide 70 that fits into the filed notch or recess 48 (Figure 5).

With reference now to Figure 8, a distal femoral resector guide 85 is mounted on the V-Frame 65 with an underside recess hole (not shown) fitting over the post 81 and with the adjustment knob 66 being rotated to adjust the stem 68 through the rotatable stem housing 69 with respect to a recess 87 in the guide 85. The guide

includes a first mount consisting of tapered laser slot 89 as well as a retractor shield 91, that locks into an "L" slot 92 and a cutting slot 93 provided to facilitate resection of the distal femur. An angular indicator 95 is also provided on the guide 85 that correlates to engraved indicia on the top surface 78. Rotation of knob 66 results in rotation of guide 85 about post 81 to adjust the relative position of guide 85 with respect to V-Frame. The indicator 95 is positioned with respect to the rotational indicia 78 on wing 75 to indicate the relative position.

With reference now to Figures 9a and 9b, it is seen that a sterile laser canister 100 includes a bracket or stem 101 that is designed to be inserted into the laser slot 89 as seen in Figure 10 with the tapered body or mount 108 under the slot 89. The stem 101 includes a swivel 103 that includes an adjustment knob 105 (Figure 10) that may be loosened to allow pivoting of the halves 104, 106 of the stem 101 about the pivot 103, whereupon the knob 105 may be tightened to lock the angular relation between the halves 104 and 106. Manipulations of the adjustment knob 66 enable one to adjust the angular position of the saw guide slot 93 as well as of the laser slot 89 with the canister 100 following along as it is mounted within the laser slot 89. Thus, the canister 100 can have its position adjusted up and down and side to side with respect to

the distal femur to facilitate aiming at the target 35 (see Figure 3).

A laser module 110 is snugly mounted within the laser canister 100 (Figure 10) and a sterile cap 111 is fastened over the opening 102 of the canister 100 to enclose the laser module therein. In the preferred embodiment of the present invention, the laser module 110 consists of a diode laser capable of emitting a highly linear beam of light, the direction of which may be adjusted in the manner explained above. Alignment is maintained between the diode laser and the canister by interengaging structures consisting of pin 112 entering into slot 113 and recesses 114 entering into pegs on the back surface of the lens face 118 module 110 (not shown). The canister 100 maintains proper alignment of the module 110 via the tapered body 108 located under the slot 89 (Figures 8 and 10). The laser module 110 has a diode laser generator 109 and an on-off switch 115 on its superior surface that turns on when it is inserted into the canister and turns off when it is removed (Figure 9b).

With the apparatus of the present invention having been described in detail, the surgical method employing this apparatus will now be described in detail.

STEP 1: THE SETUP

A) POSITION THE PATIENT ON THE OPERATIVE TABLE:

Apply the non-sterile **Operative Table Clamp 15** to the side rail of the operative table 10 as far distally as it will go in the second section of the operative table. Be sure this is on the same side as the operative knee 5. As the patient is being positioned on the operative table, use the non-sterile **Target Positioner 13** within the Operative Table Clamp 15 so that when the patient's hip is flexed to 90 degrees (Figure 1), it just touches the horizontal arm 21 of the Positioner. Move the patient up or down, or if the OR table clamp can be move cephalad without moving the patient do so. This will ensure that the Target will be correctly positioned over the patient's hip (Figure 3). Now prep the knee and drape in the usual fashion.

B) APPLY THE LASER TARGET OVER THE HIP:

Apply the sterile **Target attachment arm 18** to the operative table through a small slit through the drapes into the operative table clamp 15. Cover the slit in the drapes with a towel and tightly towel clip. Hold the **Attachment Arm 18** rigid with the adjustment holes pointing to the patient's head while the nurse tightens the Attachment Arm 18 into place into the operative table clamp. Now apply the **Horizontal Bar 27** pointing distally to the

Attachment Arm. Be sure the arm is flipped over to the correct side and tighten the adjustment screw 43 through mount 41 when the arm is about **6 inches** above the patient's abdomen (you need the height to flex the knee). Now slide the **Laser Target** system 29 on the Horizontal Bar from the open end. Be sure the center of the target portion 38 is facing distally. The actual center of the femoral head will be distal to the Hip Target leading edge by the bulk of the soft tissues on the anterior thigh.

C) ESTIMATE LOCATION OF THE FEMORAL HEAD

Adjust the **Target's** bulls-eye 36 (Figure 3) transversely 33 to be centered over the femoral head 6, approximately 2 inches medially to the anterior superior iliac spine (ASIS) (The lateral edge of the target will be directly over the ASIS). Pull the Laser Target distally 39 about 2.5 inches from the closed position. The center of the target (bulls-eye 36) will now be over the approximate center of the femoral head. It is vitally important that the leading edge of the target be pulled past the center of the femoral head by at least 2 inches to facilitate the triangulation process. If there is not enough adjustment in the target, move the operative table clamp 15 distally. This moves the entire target distally, resulting in more adjustability.

For the triangulation process to work, you will need to be able to move the operative leg at least 10 to 15 degrees (in neutral rotation) in both adduction and abduction. In the rare case of ankylosis or significant contracture of the hip, another method of alignment must be used because the center of rotation will not be obtainable by laser triangulation. Using the standard approach to total knee arthroplasty, direct your attention to the distal femur.

STEP 2: PREPARE THE DISTAL FEMUR

A) CENTRALIZE THE INTRACONDYLAR NOTCH

File the intracondylar notch with a 5/16-inch (.312) Rasp 45 (Figure 4) to form filed notch 48 (Figure 5). It is a vital importance to determine the exact center of the knee which, is usually, located slightly medial to the existing notch. Remove all the osteophytes and deepen the notch 47 with the rasp. It is important that the V-Frame 65 bottoms out on the adjacent condyles and is not held proud by the notch guide 70.

B) CREATE FLAT PLANE ON ANTERIOR FEMORAL CORTEX

File the anterior intercondylar notch 51 flush with the anterior femoral cortex with a 5/8 inch flat Rasp 49 to form flat surface 52 (Figure 5). Aim 5-7 degrees medially, which is the

direction of the femoral head. There is usually a 5-10 degree medial slope of the anterior femoral cortex so the file will take more bone off on the lateral side. The cutting guide will sit on this surface so make it as flat as possible.

5 **C) REMOVE ANTERIOR FEMORAL CONDYLES**

Use the **Anterior Condyle Resector Guide 53** (Figure 6) flush on the filed plane 52 in the intercondylar notch to resect the anterior femoral condyles 59 & 61 with the saw 57 flush with the top 55 of the guide 53. This leaves 1/8 inch of anterior condylar bone above the filed central plane 52.

10

STEP 3: APPLY HARDWARE TO DISTAL FEMUR

A) APPLY V-FRAME TO DISTAL FEMUR

Apply the **V-Frame 65** with the **Tongue 79** flush to the filed anterior femoral cortex 52 (Figure 7). Hold the 1/4 inch **Notch Rod 70** tight into the intercondylar notch recess 48 making sure the distal femoral condyles engage the **V-Frame 65** brackets 67 & 69. Pre-drill the outer cortexes only, then screw cortical screws 71 into at least two holes on each side. Make sure the **Tongue 79** is still flush with the anterior femoral cortex 52 while affixing the **V-Frame 65** to bone.

15

20

B) APPLY DISTAL FEMORAL RESECTOR TO V-FRAME

NOTE: Remove the **Tongue** 79 before proceeding with this next step. Failure to do so will impede the resection of the distal femoral condyles.

5 Apply the **Distal Femoral Resector Guide** 85 to the **V-Frame** 65 (Figure 8) and, at the same time, center the Adjustment knob stem 68 so it fits into its receptor 87 on the resector guide 85. This locks the resector guide into rough alignment with the center of the femoral head. Be sure the resector is set to the 90-degree
10 mark 95 that corresponds to the 90-degree mark of the indicia 78 on the superior surface of the V-Frame 65 (the resector indicator arm will be flush with the anterior edge of the V-Frame). Next apply the **Retractor Shield** 91 into the square slot 92 at the proximal end of the **Distal Femoral Resector Guide** 85 (Figure 8) and lock into
15 place by pushing down and distally. This retractor shield will allow good visibility of the fixation holes on the proximal arm of the Distal Femoral Resector Guide.

STEP 4: LASER SETUP & ALIGNMENT OF THE LASER

20 The Laser **Module** 110 (Figures 9a & 9b) used in this procedure cannot be sterilized! Steam sterilization would damage the delicate laser components. Therefore, the nurse inserts the non-

sterile laser module 110 into the **Sterile Laser Canister** 100. The groove 113 on the laser module must be aligned with the rib 112 on the laser canister. The laser module lines up with the lens 118 of the canister. The laser automatically illuminates when inserted, and turns off when removed by the switch 115. The surgeon closes the canister with the **sterile cap** 111. The surgeon then inserts the self-centering alignment flange 108 on the bottom of the canister into the dovetailed alignment groove 89 on the **Distal Femoral Resector Guide** 85 (Figure 10). The same Laser will later be moved to the Target mount 116 (Figure 11) for creating the Mechanical Longitudinal Axis 140 (Figure 17) to be used to check alignment throughout the case.

The **Laser Module** 110 is now locked into alignment with the **Distal Femoral Resector Guide** 85 (Figure 10) and is located at the midpoint of the knee and emits an angulated laser beam 107. What needs to be determined now is the center of the femoral head, so the Mechanical Axis of the Femur can be located accurately.

ALIGNMENT OF THE LASER

With reference to Figure 11, first apply the Laser at the Knee; line up the knee laser beam 107 to be co-linear by positioning the leg, and moving the **Target** in the transverse plane (Figure 3) until the laser beam illuminates the centerline 117

of the target 35. Be sure to hold the leg as straight as possible, parallel with the operating table. It is also very important to keep the **laser line beam** 107 positioned vertically, keeping the beam parallel with the vertical markers 42 on the **Target**. Rotation
5 of the knee will impair the accuracy of the laser beam. Therefore, it is important to always make sure the beam is parallel to the vertical lines on the front of the **Target Vertical Adjustment 38** whenever making a determination.

A word of caution before starting the triangulating process:
10 Occasionally, the laser beam will not illuminate the top of the target. This is likely in a patient with a fat thigh. In this case, use the swivel mechanism 103 on the side of the LASER to tilt the Laser beam upward until it shows on the top of the target. Also, the entire target can be moved up for a flexed thigh hitting
15 the underside of the target. Two people are needed, one to hold the leg parallel to the lines on the Target and the other to make the target adjustments.

STEP 5: LOCATING THE FEMORAL HEAD

FIRST: Determine the Longitudinal Axis of the Knee Laser

a. Line up the Knee Laser Beam 107 with the centerline 117 on the Target 35 (Figure 11).

5 This preliminary longitudinal axis will be within a few degrees of the actual longitudinal axis because of the anatomical fit of **V-Frame** 65 in the femoral intercondylar notch 48 and the tangential surfaces of the distal femoral condyles (Figure 11).

10 b. Swing the patient's flexed leg laterally (**L**) to 30 degrees and hold in place. Make a small dot on the target centerline 117 where it intersects the laser line beam 107. Label this dot **L** 120 (Figure 12).

c. Next, sweep the leg medially (**M**) to 30 degrees and hold in
15 place. Make another small dot on the target centerline 117 at the laser line beam 107 intersection. Label this dot **M** 121.
NOTE: If the points are coincident, skip ahead to Step 2 (Figure 12).

d. Move the Target longitudinally 39 so the bulls-eye 36 is on
20 the centerline of the more superior dot of the two beams, 120 for this illustration (Figure 13).

****IMPORTANT:** When dot L is superior to the M dot, the Knee laser adjustment is medial. If dot M is superior to the L dot, Knee laser adjustment is made laterally.

e. Turning to (Figure 14) Carefully measure the distance

between dot 120 L and dot 121 M with a ruler or by counting the number of 2mm lines on the built-in ruler 124 results in a distance "e" 128 divide this distance "e" by two. Draw a short midline transverse axis 130 between the dots (Figure 14).

f. Using the same midline distance in mm (one-half of "e") make a Corrective dot 131 in the direction of the correction on this midline 130 (Figure 14).

g. With the Knee Laser beam 107 returned to centerline 117 of the Target (Figures 12 and 13), dial the ADJUSTMENT KNOB 66 (Figure 15) on the **Distal Femoral Resector Guide** so the knee laser beam 107 is in a parallel position 132 with the Corrective Dot 131 (Figure 14). **This adjustment will correct the longitudinal axis of the Knee to the center of the femoral head** (Figures 14 and 16).

SECOND: Adjust the Longitudinal and Transverse Axis of the Target to Knee Laser corrective dot position:

a. Move the Target transverse axis (bulls-eye) by pulling the Target longitudinally 39 to the midpoint of the Corrective Dot 131. Measure this distance from a fixed position on the Target 133 (Figure 16).

5 b. Now make the Knee Laser beam 107 co-linear with the Target centerline 117 by moving the target transversely 33 (Figure 11). There is only one position the bulls-eye will line up with the medial and lateral radial lines 134 for this patient, which is the exact center of the femoral head
10 (Figure 16).

For the perfectionist, minor adjustments can be made here by moving the bulls-eye longitudinally and adjusting the knee laser ADJUSTMENT KNOB to match the selected radial line. After making an adjustment knob change, be sure to always make the laser co-linear
15 with the Target's centerline BY MOVING the Target transversely. If it is not lining up, simply repeat the **FIRST STEP** from your present location.

**THE BULLS-EYE OF THE LASER TARGET IS NOW POSITIONED DIRECTLY OVER THE CENTER OF THE FEMORAL HEAD, AND THE BEAM IS SHOWING THE LEG'S
20 MECHANICAL AXIS.**

Check the results by moving the Knee Laser beam 107 medially and laterally keeping it parallel with the vertical lines 42 on the front of the target and the radial lines 134 on the top of the

target. It should transverse through the center of the Target
bulls-eye 36. Do not worry if it is slightly off the center of the
bulls-eye because the Laser is extremely sensitive. If the surgeon
is within the bulls-eye **quarter inch circle**, they will be accurate
5 within a quarter of a degree of the actual center of the Femoral
Head. Spending a lot of time for perfection is not worth the
effort because it is rare to make the resection cuts that accurate.
Later minor adjustments can be made with the Mechanical Axis
Longitudinal Laser Line 140 (Figure 17) during the cementing
10 process.

The knee laser is now removed and the distal femoral cut is
made with the saw. Be sure to remove the cortical screws 71
holding the V-Frame 65 to the distal condyles before making this
resection. The V-Frame itself can also be removed if loose. The
15 laser is now moved to the second mount comprising the Target Laser
Mount 116 (Figure 17) to create the Laser Ruler 140. The laser
beam 107 will have to be tilted down via the swivel 103 by turning
locking knob 105 to shine on the entire leg creating a visual line
140 (Laser Ruler) on the leg. The remainder of the knee
20 arthroplasty is done in the conventional manner.

VALUE OF THE LONGITUDINAL AXIS RULER:

Since the hip laser is centered on the femoral head, the
straight line this laser emits is the mechanical axis for this

patient. Neutral alignment is when the laser passes through the center of the knee and the center of the ankle. Hold the leg parallel to this line to check your alignment after making all cuts with the trial components in place.

5 **HELPFUL HINTS:**

While aligning the leg on the Mechanical Axis laser-line, watch out for external rotation errors. There is a tendency for the femur to externally rotate especially in heavier patients while resting on the operative table. About 10-15 degrees of rotation
10 error can equal 1 degree of alignment error depending on the length of the patient's femur.

To avoid this error, the surgeon should stand at the end of the operating table with one hand under the distal femur rotating the femur to neutral and pushing the trial components together with
15 his abdomen. He then must align the laser beam on the center of the knee and ankle. Deviation of the laser beam from the center of the knee is the alignment error from neutral (mechanical axis).

For alignment errors less than one-degree (6-7mm), use bone cement to correct the error to neutral. For errors more than 1
20 degree (>7mm), be sure that rotation is not causing some of this error. Then check the tibial resection by re-applying the tibial resector and make sure this cut is perpendicular to the axis of the tibia. Tibial error is easily fixed by filing down the high side

evenly. Femoral errors that cannot be filed will require re-applying the distal femoral cutting jigs.

As such, an invention has been disclosed in terms of a preferred embodiment thereof including apparatus and method for its practice that fulfills each and every one of the objects of the invention as set forth hereinabove, and provides a new and useful laser triangulation of the femoral head for total knee arthroplasty alignment instruments and surgical method of use, of great novelty and utility.

Of course, various changes, modifications and alterations in the teachings of the present invention may be contemplated by those skilled in the art without departing from the intended spirit and scope thereof.

As such, it is intended that the present invention only be limited by the terms of the appended claims.